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Vákuum 2, 8, 9, 10, 11, 12, 13, 14, 18, 20
Vákuum értendő 2, 8

Z

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Electromagnetic emissions		
Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that they are used in such an environment.		
Emission tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Thopaz uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Thopaz is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 60000-3-3	Complies	

Warning – Thopaz should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Thopaz should be observed in order to verify normal operation in the configuration in which it will be used.

A Technical documentation

IEC 60601-1-2, Table 202

Electromagnetic immunity

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that they are used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 sec	< 5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Thopaz requires continued operation during power mains interruptions, it is recommended that Thopaz be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

Electromagnetic immunity

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that Thopaz is used in such an environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of Thopaz, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 2.3 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.^a</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Thopaz is used exceeds the applicable RF compliance level above, Thopaz should be observed in order to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Thopaz.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and Thopaz

Thopaz is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Thopaz can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Thopaz, as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
W	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

B Safety tests

EN IEC 62353: 2008-8 Recurrent Test and Test after Repair

For safety tests, it is assumed that Thopaz is maintained and repaired throughout its service life strictly and exclusively by Medela authorised service centres.

For Thopaz as a device in protection class II (EN IEC 60601-1), the safety tests are confined to visual inspection of the housing and mains power cord for damage. This check must be performed prior to each use.

Devices of protection class II do not have a protective earth conductor; there is therefore no need to check the earth leakage current!

The Thopaz enclosure is made entirely of insulating material. Tests of the enclosure leakage current using common measuring instruments will therefore not reveal measurable values.

The applied parts of Thopaz are connected to Thopaz via non-conductive vacuum hoses, collection jars and intermediate tubings.

Even when suctioning a conductive fluid until the overflow protection device activates (9g/L NaCl; EN ISO 10079-1;1999), measurements of the patient leakage current using common measuring instruments will not reveal measurable values.

Thopaz pumps do not have patient circuits and functional earth connections.

B Safety tests

DIN EN 62353: 2008-8

Hospital	Thopaz serial number:
	Firmware version:
Location:	Run time (Settings 2/2): _____ h

Visual checks	NO	YES
Cracks, brittle areas, discolouration on housing and display		
Covers for ports are present and can be clicked in		
Standard rail and bed holders are present and can be opened		
Seal is present, is clean and not brittle		
Release button can be pressed, canister can be clicked in		
Thopaz is clean with no dirty areas		

Measurements*	OK	Not OK
Max. vacuum (use Thopaz with 0.8 L canister and tubing set). Set vacuum to -9 kPa (-90 cmH ₂ O), refer to section X. Using a short piece of suitable tubing (7 mm ID) connect suction tube to vacuum gauge and start unit, wait 1 minute for the pump to reach equilibrium.	<input type="checkbox"/> Measured value is between -9,0 kPa (-90 cmH ₂ O) +/- 5%	<input type="checkbox"/> Measured value _____
Max. Flow (keep settings from above) Disconnect vacuum gauge and connect flow meter instead, wait 2 minutes for the pump to reach equilibrium.	<input type="checkbox"/> > 5000 mL/min	<input type="checkbox"/> Measured value _____
A warning signal (beeb and text in the display: Leak in system) may occur within the 2 minutes, otherwise keep the pump running until the warning signal (max 5 minutes).	<input type="checkbox"/> Warning 'Leak in system'	<input type="checkbox"/>
Battery Connect Thopaz to mains adapter or place it in docking station.	<input type="checkbox"/> Battery is being charged (refer to section IV)	<input type="checkbox"/> Battery is not being charged

Overall evaluation	Errors	Measures
No safety or functional technical errors were detected.	<input type="checkbox"/>	Thopaz is ready for use
Errors that are not safety relevant were detected. No direct risk.	<input type="checkbox"/> 1 _____ 2 _____ 3 _____	Thopaz is ready for use, errors must be corrected by _____
Safety relevant errors were detected.	<input type="checkbox"/> 1 _____ 2 _____ 3 _____	Do not use Thopaz until the errors are corrected. Send Thopaz to Medela Service Centre.

* Required test equipment
 Vacuum gauge ~0-10 kPa Class 2.5% or better (077.1456)
 Flow Meter ~10 L/min air 20°, 1 bar abs, Class 4% or better

Next test:

Tester / Signature:	Place / Date:
---------------------	---------------